


SimplyThick[®] Brand Pump Cleaning & Maintenance



Proper cleaning and maintenance will promote trouble-free operation of the SimplyThick[®] brand pump.


Disassembly for Cleaning

- Unscrew pump & cap assembly from container of gel.
- Pull Dip Tube ⑩ from Body ⑨.
- Separate pump by unscrewing Collar ⑤ from Body ⑨.
- Remove Spring ⑧ and Restrictor Clip ⑥.
- Separate Nozzle ② from Head ③ by pulling (it is a friction fit). Remove Cap ① from Nozzle ② by twisting and pulling (it is a friction fit). This step is optional.
- Wash with warm, soapy water. Rinse. Sanitize and dry all parts.
- Carefully reassemble pump by using diagram.

Warning: Do not clean pump in automatic dishwasher as valving could be damaged by high water temperature.

SimplyThick, LLC
200 South Hanley Road
Suite 1102
St. Louis, MO 63105

(800)205-7115 **Customer Service:** Available 24/7
(888)721-2023 **Home Office:** 8:30–4:30 CST, M–F
(800)508-2990 **Fax**
info@simplythick.com www.simplythick.com



La limpieza y el mantenimiento apropiado promoverán el funcionamiento sin problemas del dispensador de marca SimplyThick®.

Desarme para limpieza

- Desenrosque el bombeador y la tapa del dispensador de gel.
- Jale el Tubo Profundo ⑩ del Tubo Principal ⑨.
- Separe el bombeador desenroscando el Cuello ⑤ del Tubo Principal ⑨.
- Remueva el Resorte ⑧ y el Clip de Restricción ⑥.
- Separe el Pico ② de la Cabeza ③ jalándolo (se coloca a fricción). Remueva la Tapa ①, del Pico ② desenroscando y jalando (está colocado a fricción). Este paso es opcional.
- Lave con agua tibia, jabonosa. Enjuague. Desinfecte y seque todas las partes.
- Cuidadosamente vuelva a armar el dispensador usando el diagrama.

Advertencia: No lave el dispensador en el lava-vajilla automático, ya que la válvula puede dañarse con el agua a alta temperatura.

SimplyThick, LLC
200 South Hanley Road
Suite 1102
St. Louis, MO 63105

(800)205-7115 **Servicio al Cliente:**
disponible las 24 horas, los 7 días
(888)721-2023 **Oficina:** 8:30–4:30 (Zona Horaria Central)
de lunes a viernes

(800)508-2990 **Fax**
info@simplythick.com

www.simplythick.com